


Dramedilol

Chemical Properties

CAS No. : 76953-65-6
 Formula: C₂₀H₂₉N₅O₄
 Molecular Weight: 403.48
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dramedilol is a β -blocker with peripheral vasodilative activity under development by ISF (SmithKline Beecham).
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4784 mL	12.3922 mL	24.7844 mL
5 mM	0.4957 mL	2.4784 mL	4.9569 mL
10 mM	0.2478 mL	1.2392 mL	2.4784 mL
50 mM	0.0496 mL	0.2478 mL	0.4957 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481