

## Hopeachinol B

## Chemical Properties


CAS No. :

Formula: C<sub>35</sub>H<sub>24</sub>O<sub>9</sub>

Molecular Weight: 588.568

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Hopeachinol B is a useful organic compound for research related to life sciences and the catalog number is T126169.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.699 mL	8.4952 mL	16.9903 mL
5 mM	0.3398 mL	1.699 mL	3.3981 mL
10 mM	0.1699 mL	0.8495 mL	1.699 mL
50 mM	0.034 mL	0.1699 mL	0.3398 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481