


## cis-3,4-Dihydroxy-beta-ionone

## Chemical Properties

CAS No. :	875666-39-0
Formula:	C <sub>13</sub> H <sub>20</sub> O <sub>3</sub>
Molecular Weight:	224.3
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	cis-3,4-Dihydroxy-beta-ionone is a natural product of Gynura, Asteraceae. The catalog number is TN3664 and the CAS number is 875666-39-0. cis-3,4-Dihydroxy-beta-ionone can be used as a reference standard.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4583 mL	22.2916 mL	44.5831 mL
5 mM	0.8917 mL	4.4583 mL	8.9166 mL
10 mM	0.4458 mL	2.2292 mL	4.4583 mL
50 mM	0.0892 mL	0.4458 mL	0.8917 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481