


(24S)-MC 976

Chemical Properties

CAS No. :	112849-14-6
Formula:	C ₂₇ H ₄₂ O ₃
Molecular Weight:	414.62
Storage:	Powder: -20°C for 3 years Actual storage temperature shall be subject to the COA.



Biological Description

Description	(24S)-MC 976, a derivative of Vitamin D ₃ , is described as [(1R,3aR,7aR)-1-[(2S)-6,7-dihydroxy-6-methylheptan-2-yl]-3a-(2-hydroxyethyl)-7a-methyl-octahydro-1H-inden-4-yl](2E)-3-(4-methylphenyl)prop-2-enoate.
Targets(IC ₅₀)	Others,Vitamin

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4118 mL	12.0592 mL	24.1185 mL
5 mM	0.4824 mL	2.4118 mL	4.8237 mL
10 mM	0.2412 mL	1.2059 mL	2.4118 mL
50 mM	0.0482 mL	0.2412 mL	0.4824 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Strugnell S, et al. Metabolism of a cyclopropane ring-containing analog of 1 α -hydroxyvitamin D₃ in a hepatocyte cell model. Identification of 24-oxidized metabolites. *Biochemical Pharmacology* (1990), 40(2), 333-341.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481