


Tyromycin A

Chemical Properties

CAS No. :	141364-77-4
Formula:	C ₂₆ H ₃₈ O ₆
Molecular Weight:	446.58
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tyromycin A is an inhibitor of leucine and cysteine aminopeptidases.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2392 mL	11.1962 mL	22.3924 mL
5 mM	0.4478 mL	2.2392 mL	4.4785 mL
10 mM	0.2239 mL	1.1196 mL	2.2392 mL
50 mM	0.0448 mL	0.2239 mL	0.4478 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481