

## Murrayamine C

## Chemical Properties

CAS No. :	134779-19-4
Formula:	C <sub>23</sub> H <sub>25</sub> N <sub>2</sub> O <sub>2</sub>
Molecular Weight:	347.45
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Murrayamine C is a natural product for research related to life sciences. The catalog number is TN6192 and the CAS number is 73385-58-7.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8781 mL	14.3906 mL	28.7811 mL
5 mM	0.5756 mL	2.8781 mL	5.7562 mL
10 mM	0.2878 mL	1.4391 mL	2.8781 mL
50 mM	0.0576 mL	0.2878 mL	0.5756 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481