

Compound N064-0050

## Chemical Properties


CAS No. :

Formula: C<sub>40</sub>H<sub>56</sub>O<sub>4</sub>

Molecular Weight: 600.884

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound N064-0050 is a useful organic compound for research related to life sciences and the catalog number is T131199.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6642 mL	8.3211 mL	16.6423 mL
5 mM	0.3328 mL	1.6642 mL	3.3285 mL
10 mM	0.1664 mL	0.8321 mL	1.6642 mL
50 mM	0.0333 mL	0.1664 mL	0.3328 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481