

Flag Tag Immunomagnetic Beads (300 nm)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at 4°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description

TargetMol's Flag Tag Immunomagnetic Beads specifically bind to proteins with the DYKDDDDK (Flag) tag, which can be used for immunoprecipitation (IP) of proteins, protein complexes, protein-nucleic acid complexes, and other antigens. This product is suitable for antigen samples from cell lysates, cell culture supernatants, serum, ascites, etc.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481