


(RS)-3-Hydroxyphenylglycine

Chemical Properties

CAS No. :	31932-87-3
Formula:	C8H9NO3
Molecular Weight:	167.16
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(RS)-3-Hydroxyphenylglycine is a PI-linked metabotropic glutamate receptors agonist.
Targets(IC50)	Others, GluR

Solubility Information

Solubility	1eq. NaOH: 15 mg/mL (89.73 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 1 mg/mL (5.98 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.9823 mL	29.9115 mL	59.8229 mL
5 mM	1.1965 mL	5.9823 mL	11.9646 mL
10 mM	0.5982 mL	2.9911 mL	5.9823 mL
50 mM	0.1196 mL	0.5982 mL	1.1965 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Brabet I, et al. Phenylglycine derivatives discriminate between mGluR1- and mGluR5-mediated responses. Neuropharmacology. 1995 Aug;34(8):895-903.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481