

Caffeoylmalic acid

Chemical Properties


CAS No. :

Formula: C₁₃H₁₂O₈

Molecular Weight: 296.231

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Caffeoylmalic acid is a useful organic compound for research related to life sciences and the catalog number is T124908.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3758 mL	16.8788 mL	33.7576 mL
5 mM	0.6752 mL	3.3758 mL	6.7515 mL
10 mM	0.3376 mL	1.6879 mL	3.3758 mL
50 mM	0.0675 mL	0.3376 mL	0.6752 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481