


## Carpinontriol B

## Chemical Properties

CAS No. :	473451-73-9
Formula:	C <sub>19</sub> H <sub>20</sub> O <sub>6</sub>
Molecular Weight:	344.36
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Carpinontriol B has antimicrobial activity, it at 40 ug/disk caused the formation of zones of inhibition. It also has anti-inflammatory activity, it shows considerable inhibition on the production of nitric oxide and reduces the production of interleukin-6 in dose-dependent manner in RAW 264.7 cells.
Targets(IC50)	Antifection,IL Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9039 mL	14.5197 mL	29.0394 mL
5 mM	0.5808 mL	2.9039 mL	5.8079 mL
10 mM	0.2904 mL	1.452 mL	2.9039 mL
50 mM	0.0581 mL	0.2904 mL	0.5808 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481