


2',4'-Dihydroxy-3'-methoxyacetophenone

Chemical Properties

CAS No. :	62615-26-3
Formula:	C ₉ H ₁₀ O ₄
Molecular Weight:	182.17
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2',4'-Dihydroxy-3'-methoxyacetophenone is a natural product for research related to life sciences. The catalog number is TN2717 and the CAS number is 62615-26-3.
Targets(IC50)	ADC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4894 mL	27.4469 mL	54.8938 mL
5 mM	1.0979 mL	5.4894 mL	10.9788 mL
10 mM	0.5489 mL	2.7447 mL	5.4894 mL
50 mM	0.1098 mL	0.5489 mL	1.0979 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481