

BDP 576/589 tetrazine


Chemical Properties

CAS No. :

Formula: C₂₆H₂₃BF₂N₈O

Molecular Weight: 512.32

Storage: Keep away from direct sunlight
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	BDP 576/589 tetrazine shares similar characteristics with other dyes in the BDP family, distinguished by its terminal tetrazine group.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9519 mL	9.7595 mL	19.5191 mL
5 mM	0.3904 mL	1.9519 mL	3.9038 mL
10 mM	0.1952 mL	0.976 mL	1.9519 mL
50 mM	0.039 mL	0.1952 mL	0.3904 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481