

BDP 558/568 DBCO

Chemical Properties

CAS No. :


Formula: C₃₇H₃₃BF₂N₄O₂S

Molecular Weight: 646.56

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	BDP 558/568 DBCO is a yellow dye featuring a terminal DBCO group.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5466 mL	7.7332 mL	15.4665 mL
5 mM	0.3093 mL	1.5466 mL	3.0933 mL
10 mM	0.1547 mL	0.7733 mL	1.5466 mL
50 mM	0.0309 mL	0.1547 mL	0.3093 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481