


## Idazoxan

## Chemical Properties

CAS No. :	79944-58-4
Formula:	C <sub>11</sub> H <sub>12</sub> N <sub>2</sub> O <sub>2</sub>
Molecular Weight:	204.23
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Idazoxan is an alpha2 receptor antagonist that also exhibits activity at imidazoline [I1] and [I2] receptors, playing a role in modulating dopamine release.
Targets(IC50)	Adrenergic Receptor, Imidazoline Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8964 mL	24.4822 mL	48.9644 mL
5 mM	0.9793 mL	4.8964 mL	9.7929 mL
10 mM	0.4896 mL	2.4482 mL	4.8964 mL
50 mM	0.0979 mL	0.4896 mL	0.9793 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481