

GJN75561

Chemical Properties

CAS No. : 2253775-56-1


Formula: C₂₁H₂₆D₆ 18FNO₃

Molecular Weight: 370.52

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

GJN75561 is a VMAT2 imaging agent.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6989 mL	13.4945 mL	26.9891 mL
5 mM	0.5398 mL	2.6989 mL	5.3978 mL
10 mM	0.2699 mL	1.3495 mL	2.6989 mL
50 mM	0.054 mL	0.2699 mL	0.5398 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481