


Gestacalone

Chemical Properties

CAS No. : 19291-69-1
 Formula: C₂₃H₂₇ClO₂
 Molecular Weight: 370.91
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Gestacalone is a progesterone receptor agonist, designed specifically for the activation of these receptors.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6961 mL	13.4804 mL	26.9607 mL
5 mM	0.5392 mL	2.6961 mL	5.3921 mL
10 mM	0.2696 mL	1.348 mL	2.6961 mL
50 mM	0.0539 mL	0.2696 mL	0.5392 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481