

RV6153

Chemical Properties


CAS No. : 1628139-03-6

Formula: C₃₆H₃₁F₄N₇O₅

Molecular Weight: 717.67

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	RV6153 is a PI3K inhibitor.
-------------	-----------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3934 mL	6.967 mL	13.934 mL
5 mM	0.2787 mL	1.3934 mL	2.7868 mL
10 mM	0.1393 mL	0.6967 mL	1.3934 mL
50 mM	0.0279 mL	0.1393 mL	0.2787 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481