


Tolgabide

Chemical Properties

CAS No. :	86914-11-6
Formula:	C ₁₈ H ₁₈ Cl ₂ N ₂ O ₂
Molecular Weight:	365.25
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tolgabide is a gabamimetic compound with properties that include antiepileptic, anticonvulsant, and anti-dyskinetic activities.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7379 mL	13.6893 mL	27.3785 mL
5 mM	0.5476 mL	2.7379 mL	5.4757 mL
10 mM	0.2738 mL	1.3689 mL	2.7379 mL
50 mM	0.0548 mL	0.2738 mL	0.5476 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481