



SDZ 210-086 maleate

Chemical Properties

CAS No. : 104925-23-7
 Formula: C₁₃H₂₁NO₆
 Molecular Weight: 287.31
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	SDZ 210-086, an oral muscarinic agonist, has been shown to reduce rapid eye movement (REM) latency, increase the percentage of REM sleep, and extend the total duration of REM sleep. Additionally, SDZ 210-086 significantly decreases the time spent in slow-wave sleep.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4806 mL	17.4028 mL	34.8056 mL
5 mM	0.6961 mL	3.4806 mL	6.9611 mL
10 mM	0.3481 mL	1.7403 mL	3.4806 mL
50 mM	0.0696 mL	0.3481 mL	0.6961 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481