


Quinotolast

Chemical Properties

CAS No. :	101193-40-2
Formula:	C17H12N6O3
Molecular Weight:	348.32
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Quinotolast (FK021) is an orally active mast cell stabilizer with cytoprotective properties on the gastric mucosa. It exhibits both antiallergic and antiulcer functions.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8709 mL	14.3546 mL	28.7092 mL
5 mM	0.5742 mL	2.8709 mL	5.7418 mL
10 mM	0.2871 mL	1.4355 mL	2.8709 mL
50 mM	0.0574 mL	0.2871 mL	0.5742 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481