


OUL243

Chemical Properties

CAS No. : 105873-33-4
Formula: C₉H₈N₄S
Molecular Weight: 204.25
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description OUL243 is a potent and selective inhibitor of mono-ADP-ribosyltransferase (ARTs), effectively suppressing PARP7, PARP10, PARP11, PARP12, PARP14, and PARP15.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.896 mL	24.4798 mL	48.9596 mL
5 mM	0.9792 mL	4.896 mL	9.7919 mL
10 mM	0.4896 mL	2.448 mL	4.896 mL
50 mM	0.0979 mL	0.4896 mL	0.9792 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481