


Preimplantation factor

Chemical Properties

CAS No. :	485818-40-4
Formula:	C67H115N21O23S
Molecular Weight:	1614.82
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Preimplantation Factor (PIF-1) exhibits potential activity against graft-versus-host disease (GVHD).
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6193 mL	3.0963 mL	6.1926 mL
5 mM	0.1239 mL	0.6193 mL	1.2385 mL
10 mM	0.0619 mL	0.3096 mL	0.6193 mL
50 mM	0.0124 mL	0.0619 mL	0.1239 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481