


Prosopine

Chemical Properties

CAS No. : 14058-38-9
 Formula: C₁₈H₃₇N₃O₃
 Molecular Weight: 315.49
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	PROSOPINE is an alkaloid derived from a leguminous shrub.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1697 mL	15.8484 mL	31.6967 mL
5 mM	0.6339 mL	3.1697 mL	6.3393 mL
10 mM	0.317 mL	1.5848 mL	3.1697 mL
50 mM	0.0634 mL	0.317 mL	0.6339 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481