


Prizidilol hydrochloride

Chemical Properties

CAS No. : 72849-00-4
 Formula: C₁₇H₂₉Cl₂N₅O₃
 Molecular Weight: 422.35
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Prizidilol (SKF 92657) is an aryl pyrazoline hydrazine derivative known for its antihypertensive effects.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3677 mL	11.8385 mL	23.677 mL
5 mM	0.4735 mL	2.3677 mL	4.7354 mL
10 mM	0.2368 mL	1.1839 mL	2.3677 mL
50 mM	0.0474 mL	0.2368 mL	0.4735 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481