


Flotegatide

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 1010702-99-4 |
| Formula: | C41H60FN13O13 |
| Molecular Weight: | 961.99 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|-------------|---|
| Description | FLOTEGATIDE is a cyclic tripeptide containing an arginine-glycine-aspartic acid (R-G-D) sequence. |
|-------------|---|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|------------|
| 1 mM | 1.0395 mL | 5.1976 mL | 10.3951 mL |
| 5 mM | 0.2079 mL | 1.0395 mL | 2.079 mL |
| 10 mM | 0.104 mL | 0.5198 mL | 1.0395 mL |
| 50 mM | 0.0208 mL | 0.104 mL | 0.2079 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481