


## Loxanast

## Chemical Properties

CAS No. : 69915-62-4  
 Formula: C<sub>14</sub>H<sub>26</sub>O<sub>2</sub>  
 Molecular Weight: 226.36  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Loxanast is an inhibitor employed in animal models to suppress delayed-type hypersensitivity reactions. It effectively inhibits the chemotactic activity of macrophages.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4177 mL	22.0887 mL	44.1774 mL
5 mM	0.8835 mL	4.4177 mL	8.8355 mL
10 mM	0.4418 mL	2.2089 mL	4.4177 mL
50 mM	0.0884 mL	0.4418 mL	0.8835 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481