


## Isomolpan

## Chemical Properties

CAS No. :	107320-86-5
Formula:	C <sub>15</sub> H <sub>21</sub> NO <sub>2</sub>
Molecular Weight:	247.33
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	ISOMOLPAN is a selective agonist for D3/D2 receptors.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0432 mL	20.2159 mL	40.4318 mL
5 mM	0.8086 mL	4.0432 mL	8.0864 mL
10 mM	0.4043 mL	2.0216 mL	4.0432 mL
50 mM	0.0809 mL	0.4043 mL	0.8086 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481