


L-731988

## Chemical Properties

CAS No. : 251922-77-7  
Formula: C<sub>15</sub>H<sub>12</sub>FNO<sub>4</sub>  
Molecular Weight: 289.26  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	L-731988 functions as an anti-HIV agent by inhibiting strand transfer, blocking integration, and hindering the replication of HIV-1 within cells.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4571 mL	17.2855 mL	34.571 mL
5 mM	0.6914 mL	3.4571 mL	6.9142 mL
10 mM	0.3457 mL	1.7285 mL	3.4571 mL
50 mM	0.0691 mL	0.3457 mL	0.6914 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481