


IGGi-11

Chemical Properties

CAS No. : 303059-39-4
 Formula: C₂₃H₂₈N₂O₈S₂
 Molecular Weight: 524.61
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	IGGi-11, an inhibitor of heterotrimeric G protein signaling, atypically activates pathways by disrupting the interaction between G α and GIV/Girdin, thereby reducing or preventing the metastasis of tumor cells.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9062 mL	9.5309 mL	19.0618 mL
5 mM	0.3812 mL	1.9062 mL	3.8124 mL
10 mM	0.1906 mL	0.9531 mL	1.9062 mL
50 mM	0.0381 mL	0.1906 mL	0.3812 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481