


## Mukonidine

## Chemical Properties

CAS No. :	24949-00-6
Formula:	C <sub>14</sub> H <sub>11</sub> NO <sub>3</sub>
Molecular Weight:	241.24
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Mukonidine is a naturally occurring carbazole alkaloid.
Targets(IC50)	Antibacterial

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1452 mL	20.7262 mL	41.4525 mL
5 mM	0.829 mL	4.1452 mL	8.2905 mL
10 mM	0.4145 mL	2.0726 mL	4.1452 mL
50 mM	0.0829 mL	0.4145 mL	0.829 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481