


Phencarbamide napadisilate

Chemical Properties

CAS No. : 72017-58-4
 Formula: C₂₉H₃₂N₂O₇S₃
 Molecular Weight: 616.77
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Phencarbamide napadisilate is an antispasmodic and anticholinergic compound used during labor.
-------------	------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6213 mL	8.1067 mL	16.2135 mL
5 mM	0.3243 mL	1.6213 mL	3.2427 mL
10 mM	0.1621 mL	0.8107 mL	1.6213 mL
50 mM	0.0324 mL	0.1621 mL	0.3243 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481