


## Panepophenanthrin

## Chemical Properties

CAS No. :	409064-57-9
Formula:	C <sub>22</sub> H <sub>28</sub> O <sub>8</sub>
Molecular Weight:	420.45
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PANÉPOPHENANTHRIN is a metabolite of Panapophenanthrin, derived from a novel oligocyclic diterpenoid isolated from <i>Panus strigellus</i> . PANÉPOPHENANTHRIN has the ability to inhibit ubiquitin-activating enzyme.
Targets(IC50)	E1/E2/E3 Enzyme

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3784 mL	11.892 mL	23.784 mL
5 mM	0.4757 mL	2.3784 mL	4.7568 mL
10 mM	0.2378 mL	1.1892 mL	2.3784 mL
50 mM	0.0476 mL	0.2378 mL	0.4757 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481