


Prodipine

Chemical Properties

CAS No. :	31314-38-2
Formula:	C ₂₀ H ₂₅ N
Molecular Weight:	279.42
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Prodipine (also known as BY-101) is a diphenylpiperidine derivative used as a compound for treating Parkinson's disease.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5788 mL	17.8942 mL	35.7884 mL
5 mM	0.7158 mL	3.5788 mL	7.1577 mL
10 mM	0.3579 mL	1.7894 mL	3.5788 mL
50 mM	0.0716 mL	0.3579 mL	0.7158 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481