


Rilozarone

Chemical Properties

CAS No. :	79282-39-6
Formula:	C ₃₂ H ₃₆ BrClN ₂ O ₂
Molecular Weight:	596
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Rilozarone has been evaluated as a coronary vasodilator and antiarrhythmic agent.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6779 mL	8.3893 mL	16.7785 mL
5 mM	0.3356 mL	1.6779 mL	3.3557 mL
10 mM	0.1678 mL	0.8389 mL	1.6779 mL
50 mM	0.0336 mL	0.1678 mL	0.3356 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481