


Sodium picofosfate

Chemical Properties

CAS No. :	36175-05-0
Formula:	C ₁₈ H ₁₃ NNa ₄ O ₈ P ₂
Molecular Weight:	525.2
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sodium picofosfate is a contact irritant laxative primarily used to alleviate constipation or prepare the colon for a colonoscopy or surgery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.904 mL	9.5202 mL	19.0404 mL
5 mM	0.3808 mL	1.904 mL	3.8081 mL
10 mM	0.1904 mL	0.952 mL	1.904 mL
50 mM	0.0381 mL	0.1904 mL	0.3808 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481