


Quincarbate

Chemical Properties

CAS No. : 54340-59-9
Formula: C₁₇H₁₈ClNO₆
Molecular Weight: 367.78
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Quincarbate is a quinoline derivative with diuretic properties.
-------------	-----------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.719 mL	13.5951 mL	27.1902 mL
5 mM	0.5438 mL	2.719 mL	5.438 mL
10 mM	0.2719 mL	1.3595 mL	2.719 mL
50 mM	0.0544 mL	0.2719 mL	0.5438 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481