


## Quinazosin

## Chemical Properties

CAS No. :	15793-38-1
Formula:	C17H23N5O2
Molecular Weight:	329.4
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Quinazosin is an antihypertensive agent.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0358 mL	15.1791 mL	30.3582 mL
5 mM	0.6072 mL	3.0358 mL	6.0716 mL
10 mM	0.3036 mL	1.5179 mL	3.0358 mL
50 mM	0.0607 mL	0.3036 mL	0.6072 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481