


Sulfachrysoidine sodium

Chemical Properties

CAS No. :	857537-62-3
Formula:	C ₁₃ H ₁₂ N ₅ NaO ₄ S
Molecular Weight:	357.32
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Sulfachrysoidine, a lipophilic azo dye derived from sulfamylchrysoidine, exhibits antibacterial properties.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7986 mL	13.9931 mL	27.9861 mL
5 mM	0.5597 mL	2.7986 mL	5.5972 mL
10 mM	0.2799 mL	1.3993 mL	2.7986 mL
50 mM	0.056 mL	0.2799 mL	0.5597 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481