


## Tocophersolan

## Chemical Properties

CAS No. :	30999-06-5
Formula:	C <sub>35</sub> H <sub>58</sub> O <sub>6</sub>
Molecular Weight:	574.83
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Tocophersolan, a synthetic polymer and form of vitamin E, can independently form micelles because of its amphiphilic nature, allowing for easy absorption by intestinal epithelial cells even in the absence of bile salts.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7396 mL	8.6982 mL	17.3964 mL
5 mM	0.3479 mL	1.7396 mL	3.4793 mL
10 mM	0.174 mL	0.8698 mL	1.7396 mL
50 mM	0.0348 mL	0.174 mL	0.3479 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481