


## S-Tenatoprazole sodium

## Chemical Properties

CAS No. :	705968-89-4
Formula:	C <sub>16</sub> H <sub>17</sub> N <sub>4</sub> NaO <sub>3</sub> S
Molecular Weight:	368.39
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	S-Tenatoprazole sodium: A promising new proton pump inhibitor with favorable safety pharmacology characteristics.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7145 mL	13.5726 mL	27.1451 mL
5 mM	0.5429 mL	2.7145 mL	5.429 mL
10 mM	0.2715 mL	1.3573 mL	2.7145 mL
50 mM	0.0543 mL	0.2715 mL	0.5429 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481