


Tazadolene

Chemical Properties

CAS No. : 87936-75-2
 Formula: C₁₆H₂₁N
 Molecular Weight: 227.35
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tazadolene is a non-opioid analgesic with antidepressant properties.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3985 mL	21.9925 mL	43.985 mL
5 mM	0.8797 mL	4.3985 mL	8.797 mL
10 mM	0.4399 mL	2.1993 mL	4.3985 mL
50 mM	0.088 mL	0.4399 mL	0.8797 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481