


Tametraline HCl

Chemical Properties

CAS No. : 52760-47-1
 Formula: C₁₇H₂₀ClN
 Molecular Weight: 273.8
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tametraline (also known as CP-24,441) is a norepinephrine-dopamine reuptake inhibitor developed for the treatment of depression.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6523 mL	18.2615 mL	36.523 mL
5 mM	0.7305 mL	3.6523 mL	7.3046 mL
10 mM	0.3652 mL	1.8262 mL	3.6523 mL
50 mM	0.073 mL	0.3652 mL	0.7305 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481