


## Suvadronabinol

## Chemical Properties

CAS No. : 1225194-84-2  
 Formula: C<sub>30</sub>H<sub>43</sub>NO<sub>6</sub>  
 Molecular Weight: 513.67  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	Suvadronabinol is a cannabinoid receptor agonist.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9468 mL	9.7339 mL	19.4678 mL
5 mM	0.3894 mL	1.9468 mL	3.8936 mL
10 mM	0.1947 mL	0.9734 mL	1.9468 mL
50 mM	0.0389 mL	0.1947 mL	0.3894 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481