


Tipredane

Chemical Properties

CAS No. : 85197-77-9
 Formula: C₂₂H₃₁F₀S₂
 Molecular Weight: 410.61
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tipredane is a sulfur-containing steroid compound developed for the topical treatment of human skin disorders.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4354 mL	12.177 mL	24.354 mL
5 mM	0.4871 mL	2.4354 mL	4.8708 mL
10 mM	0.2435 mL	1.2177 mL	2.4354 mL
50 mM	0.0487 mL	0.2435 mL	0.4871 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481