


Teglicar HCl

Chemical Properties

CAS No. :	908566-80-3
Formula:	C ₂₂ H ₄₆ ClN ₃ O ₃
Molecular Weight:	436.07
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Teglicar (also known as ST1326) is a selective reversible inhibitor of carnitine palmitoyltransferase 1 (L-CPT1), the liver isoform. It reduces ketogenesis and glucose production, decreases gluconeogenesis, and enhances glucose homeostasis.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2932 mL	11.466 mL	22.9321 mL
5 mM	0.4586 mL	2.2932 mL	4.5864 mL
10 mM	0.2293 mL	1.1466 mL	2.2932 mL
50 mM	0.0459 mL	0.2293 mL	0.4586 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481