


Theophylline ephedrine

Chemical Properties

CAS No. :	15766-94-6
Formula:	C17H23N5O3
Molecular Weight:	345.4
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Theophylline ephedrine is a pharmacologically active salt used for the treatment of obstructive ventilatory disorders.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8952 mL	14.476 mL	28.9519 mL
5 mM	0.579 mL	2.8952 mL	5.7904 mL
10 mM	0.2895 mL	1.4476 mL	2.8952 mL
50 mM	0.0579 mL	0.2895 mL	0.579 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481