


Tolprocarb

Chemical Properties

CAS No. :	911499-62-2
Formula:	C ₁₆ H ₂₁ F ₃ N ₂ O ₃
Molecular Weight:	346.35
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	TOLPROCARB is a fungicide effective against oomycetes, functioning by inhibiting polyketide synthase in the melanin biosynthesis pathway.
Targets(IC50)	Antifungal

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8873 mL	14.4363 mL	28.8725 mL
5 mM	0.5775 mL	2.8873 mL	5.7745 mL
10 mM	0.2887 mL	1.4436 mL	2.8873 mL
50 mM	0.0577 mL	0.2887 mL	0.5775 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481