


YEATS Inhibitor 7

Chemical Properties

CAS No. :	2803461-76-7
Formula:	C17H19N5O
Molecular Weight:	309.37
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	YEATS Inhibitor 7 is an inhibitor targeting the ENL YEATS domain, used in the treatment of leukemia.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2324 mL	16.1619 mL	32.3238 mL
5 mM	0.6465 mL	3.2324 mL	6.4648 mL
10 mM	0.3232 mL	1.6162 mL	3.2324 mL
50 mM	0.0646 mL	0.3232 mL	0.6465 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481