


Trimedoxime dichloride

Chemical Properties

CAS No. :	3613-82-9
Formula:	C ₁₅ H ₁₈ Cl ₂ N ₄ O ₂
Molecular Weight:	357.24
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Trimedoxime is a cholinesterase activator frequently utilized as an antidote for alkyl phosphate poisoning.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7992 mL	13.9962 mL	27.9924 mL
5 mM	0.5598 mL	2.7992 mL	5.5985 mL
10 mM	0.2799 mL	1.3996 mL	2.7992 mL
50 mM	0.056 mL	0.2799 mL	0.5598 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481